

ISO1000R RECUMBENT BIKE

Key Features & Benefits

- Comfortable, oversized seat featuring exceptional lumbar support
- Oversized, self-righting pedals with easily adjustable straps
- Easy step-through access to allow entry without having to lift the leg up
- Heavy-duty, 3-piece proprietary crank system with forward-only resistance
- Only 6 watts of resistance needed to start and maintain programme function, ideal in rehab applications or for the 'not yet fit' user
- Certified as a Medical Device



Specification	ISO1000R
Workout Position	Recumbent
Console	Standard
Standard Seat	Tall Back Fixed
Pedals	Oversized • adjustable footstraps
Cranks	Fixed cranks as standard • adjustable available as an option
Medical CE Ila Certification	Yes
Bi-directional Resistance	Yes
Resistance	200 levels • adjustable in 0.1 increments
Programmes	Manual • Iso-Strength • Heart Rate • Constant Work • Random • Hill Profiles • Fit-Quik • Power Burst™ • Heart-Fit Test • Heart-Fit Training • Power-Fit Test • Power-Fit Training
Feedback	Heart Rate • Time • RPM • Calories • Distance • Level • Watts • METS
Power	Self-generating from just 6 watts with auto recharge and battery back up
Length	150cm
Width	66cm
Height	132cm
Weight	95kg
Max. User Capacity	204kg
Colour	Silver powder coated frame and grey metallic covers with UV protection
Warranty	3 years parts



Medical CE Ila certification and BS EN ISO13485:2008 for use in hospitals and clinics