

ISO7000 UPRIGHT BIKE

Key Features & Benefits

- Large user seat with forward/backwards and vertical adjustments for all user heights
- Oversized, self-righting pedals with easily adjustable straps
- Easy step-through access to allow entry without lifting the leg up and over central console
- Heavy-duty, 3-piece proprietary crank system with bi-directional resistance
- Only 6 watts of resistance needed to start and maintain programme function, ideal in rehab applications or for the 'not yet fit' user
- Certified as a Medical Device



Specification	ISO7000
Workout Position	Upright
Console	Standard
Standard Seat	Vertical & Horizontal Adjustments
Pedals	Oversized • self-righting
Cranks	Fixed cranks as standard • adjustable available as an option
Medical CE Ila Certification	Yes
Bi-directional Resistance	Yes
Resistance	200 levels • adjustable in 0.1 increments
Programmes	Manual • Iso-Strength • Heart Rate • Constant Work • Random • Hill Profiles • Fit-Quik • Power Burst™ • Heart-Fit Test • Heart-Fit Training • Power-Fit Test • Power-Fit Training
Feedback	Heart Rate • Time • RPM • Calories • Distance • Level • Watts • METS
Power	Self-generating from just 6 watts with auto recharge and battery back up
Length	147cm
Width	61cm
Height	134cm
Weight	51kg
Max. User Capacity	159kg
Colour	Silver powder coated frame and grey metallic covers with UV protection
Warranty	3 years parts



Medical CE Ila certification and BS EN ISO13485:2008 for use in hospitals and clinics