

## SXT7000e<sup>2</sup> EASY-ENTRY ELLIPTICAL

### Key Features & Benefits

- Easy-entry package includes sturdy rear step to reduce 'step-up' height by half, allowing access to more users and the full side handrails provide safe and easy entry on to the elliptical
- Comfortable stride length ensures natural total body movement, proper posture, optimised torso rotation and core recruitment
- Orthopedic footbeds with front/side safety edges and patented Bio-Flex™ technology provide better circulation and eliminate hot spots and numbness common with other ellipticals
- Bi-directional resistance helps create balance between reciprocal muscle groups
- The pedal movement simulates walking on an uneven surface to help recruit stabilising muscles
- Multi-position hand grips for extra comfort and elliptical has low centre of gravity for increased user safety and stability
- Starting wattage of 6 with low minimum RPM to begin and maintain programme function, ideal in rehab applications or for the 'not yet fit' user
- Certified as a Medical Device



| Specification                | SXT7000e <sup>2</sup>   |
|------------------------------|---|
| Workout Position             | Upright   |
| Console                      | Standard  |
| Standard Seat                | n/a   |
| Hangrips & Footbeds          | Dual-position handgrips and orthopedic footbeds with patented Bio-Flex™ technology  |
| Easy Entry Step & Handrails  | Yes   |
| Medical CE Ila Certification | Yes   |
| Contact Heart Rate           | Yes   |
| Resistance                   | 200 levels of bi-directional resistance adjustable in 0.1 increments  |
| Programmes                   | Manual • Iso-Strength • Heart Rate • Constant Work • Random • Hill Profiles • Fit-Quik • Power Burst™ • Heart-Fit Test • Heart-Fit Training • Power-Fit Test • Power-Fit Training |
| Feedback                     | Heart Rate • Time • RPM • Calories • Distance • Level • Watts • METS  |
| Power                        | Self-generating from just 6 watts with auto recharge and battery back up  |
| Length                       | 190cm   |
| Width                        | 76cm  |
| Height                       | 160cm   |
| Weight                       | 138kg   |
| Max. User Capacity           | 193kg   |
| Colour                       | Silver powder coated frame and grey metallic covers with UV protection  |
| Warranty                     | 3 years parts   |



Medical CE Ila certification and BS EN ISO13485:2008 for use in hospitals and clinics