

PRO2[®] TOTAL BODY EXERCISER

Key Features & Benefits

- Isolate training to be upper body only, lower body only or a total body workout
- Adjustable console angle to accommodate users at all positions
- Adjustable and reversible arm cranks for a greater range of movement and exercise variety
- Vertical seat adjustment with gas assisted lift and swivel
- Seat is removable for use with a wheelchair
- Only 6 watts of resistance needed to start and maintain programme function, ideal in rehab applications or for the 'not yet fit' user
- Certified as a Medical Device



Specification	PRO2 [®]
Workout	Total Body
Head	Fixed
Cranks	Adjustable
Console	Standard
Standard Seat	Swivel Seat
Removable Seat	Yes
Wheelchair Platform	Optional
Medical CE IIa Certification	Yes
Resistance	200 levels of bi-directional resistance adjustable in 0.1 increments
Programmes	Manual • Iso-Strength • Heart Rate • Constant Work • Random • Hill Profiles • Fit-Quik • Power Burst™ • Heart-Fit Test • Heart-Fit Training • Power-Fit Test • Power-Fit Training
Feedback	Heart Rate • Time • RPM • Calories • Distance • Level • Watts • METS
Power	Self-generating from just 6 watts with auto recharge and battery back up
Length	155cm
Width	76cm
Height	156cm
Weight	100kg
Max. User Capacity	204kg
Colour	Silver powder coated frame and grey metallic covers with UV protection
Warranty	3 years parts



Medical CE IIa certification and BS EN ISO13485:2008 for use in hospitals and clinics