

## AC5000 TREADMILL

### Key Features & Benefits

- Slow starting speed of just 0.5kph up to 19kph in 0.1 increments for extra safety and for those with mobility limitations
- Electronic elevation from 0-15% in 0.5 increments of adjustment
- Oversized belt with large user capacity and a comfortable walking/running surface
- Reversible, eucalyptus deck is very low maintenance - no lubrication required
- Side handrail switches for quick and easy speed and elevation changes
- Commercial grade 3" crowned rollers and powerful 3 HP AC motor with inverter
- Dual emergency stop button for added safety and integrated ready rack and water bottle holder
- Optional full-length side handrails and treadmill step are available
- Certified as a Medical Device



Specification	AC5000
Platform & Belt	Extra wide with large low static belt (56cm x 157cm)
Deck	Eucalyptus - maintenance-free no lubrication
Side Switches	Yes
Speed Increments	0.5kph
Speed Range	0.5-19kph
Reverse Speed	No
Elevation Range	0% to 15%
Medical CE Ila Certification	Yes
Full Side Handrails + Step	Optional
Programmes	Manual • Heart Rate • Programme Profiles
Feedback	Speed • Elevation • Heart Rate • Time • Calories Distance • Pace • METS
Power	3HP Motor with inverter, mains power required
Length	207cm
Width	80cm
Height	160cm
Weight	220kg
Max. User Capacity	227kg
Colour	Silver powder coated frame and black/silver metallic covers with UV protection
Warranty	3 years parts



Medical CE Ila certification and BS EN ISO13485:2008 for use in hospitals and clinics