

AC5000M TREADMILL

Key Features & Benefits

- Very slow starting speed from 0.2kph up to 16kph in 0.1 increments for extra safety and for those with mobility limitations and belt can also travel in the reverse direction up to speeds of 6kph for use in rehab applications
- Electronic elevation from -3% to 12% in 0.5 increments of adjustment - walking on the decline is often required in rehab applications
- Oversized belt with large user capacity and a comfortable walking/running surface
- Reversible, eucalyptus deck is very low maintenance - no lubrication required
- Side handrail switches for quick and easy speed and elevation changes
- Commercial grade 3" crowned rollers and powerful 3 HP AC motor with inverter
- Dual emergency stop button for added safety
- Optional full-length side handrails (pictured) and treadmill step are also available



Specification	AC5000M
Platform & Belt	Extra wide with large low static belt (56cm x 157cm)
Deck	Eucalyptus - maintenance-free no lubrication
Side Switches	Yes
Speed Increments	0.2kph
Speed Range	0.2-16kph
Reverse Speed	up to 6kph
Elevation Range	-3% to 12%
Medical CE IIa Certification	No
Full Side Handrails + Step	Optional
Programmes	Manual • Heart Rate • Programme Profiles
Feedback	Speed • Elevation • Heart Rate • Time • Calories Distance • Pace • METS
Power	3HP Motor with inverter, mains power required
Length	207cm
Width	80cm
Height	160cm
Weight	220kg
Max. User Capacity	227kg
Colour	Silver powder coated frame and black/silver metallic covers with UV protection
Warranty	3 years parts